



## MANUAL BOOK

# UV ROOM STERILIZER MUV-01



## ✓ INTRODUCTION

One of the unit used for the room sterilization using ultraviolet lamp. The ultraviolet lamp can be use in the operating room, ICU and others that require a room sterile conditions with the room size up to 5 x 5 m<sup>2</sup> and ultraviolet lamp life up to 3000 hours. For supporting the conditions the factory planned to create a lamp that can meet the requirement. We are very concerned about the quality of our goods on the market in order to be well received in any hospital or clinic.

For that we create a tool that we called UV ROOM STERILIZER with type MUV-01.



Gambar MUV-01

## ✓ STANDARD OPERASIONAL PROCEDURE

The unit is very easy to operate, because its shape is simple but still elegant, here's how the operation of the appliance:

### A. MANUAL PROCEDURE

1. Install the unit correctly
2. Connect the power cord into the power supply
3. Select the manual operation with manual selection button
4. Press POWER
5. Select the light you want to run with the push of a button 1, 2, 3, 4, 5 or 1-5
6. Set a timer to set the duration of sterilization
7. Once the time is up, it indicates the sterilization process has been completed and the lights to be extinguished

### B. REMOTE PROCEDURE

1. Install the unit correctly
2. Connect the power cord into the power supply
3. Select the remote operation with remote selection button
4. Press POWER
5. Select the light you want to run with the push of a button 1, 2, 3, 4, 5 or 1-5
6. Set a timer to set the duration of sterilization
7. Press the A / B / C to turn on the lights
8. Once the time is up, it indicates the sterilization process has been completed and the lights be extinguished. When suddenly want to turn off the lights, press the D button on the remote

## ✓ MAINTENANCE AND REPAIR

### Repair

No	Damage	Possible Cause	Action
1	The lamp may not be lit	No power supply Lamp broke Ballast broke Timer broke	Cek power supply Replace the lamp Replace ballast Replace timer
2	Remote tidak berfungsi	Broken remote key	Fix the button

### Maintenance

Clean the body with a soft cloth, and when using a liquid, use a gentle liquid, like soap.

File No.: SPA-BM/PROD-11

File Ver.: 0.1

Release Date: 2016.05.01 rev00