


**Statement**  
This manual will help you understand the operation and maintenance of the product properly. Please note that the product must be used correctly and must comply with this manual. Users who do not comply with this manual may result in the unit not functioning properly or an accident in which PT. Sinko Prima Alloy Instrument. (hereinafter PT. Sinko Prima Alloy) cannot be held responsible. PT. Sinko Prima Alloy owns the copyright of this manual. Without prior written approval from PT. Sinko Prima Alloy. All materials related to this manual may not be photocopied, reproduced, or translated into other languages.  
Materials are protected by copyright law, but not limited to confidential information such as technical information and patent information contained in this manual, users may not provide information to third parties that are irrelevant.  
The user should understand that nothing in this manual provides express or implicit permission to the user. All rights or licenses to the intellectual property of PT. Sinko Prima Alloy.  
PT. Sinko Prima Alloy reserves the right to modify, update and explain this manual.

**Responsibility of the Manufacturer**  
PT. Sinko Prima Alloy only considers itself responsible for any effect on the safety, constraints and equipment performance of the unit, if:  
The assembly, extension, re-adjustment, modification or repair operations are carried out by a person authorized by PT. Sinko Prima Alloy and electrical installations in accordance with National Standards. The instrument is used according to the instructions for use. PT. Sinko Prima Alloy will provide requests for circuit diagrams, part lists, descriptions, calibration instructions or other information on individual repairs to repair components of the unit appointed by PT. Sinko Prima Alloy as a component that can be repaired by service personnel.  
**Product Information**  
**Product Name:** Ultrasonic Pocket Doppler  
**Model:** DP1  
**Terms in the Manual**  
This manual is intended to provide a basic concept of security measures.

**Warning**  
**A Warning** label suggests a specific action or situation that results in personal injury or death.  
**Attention**  
**A Caution** label suggests an action or situation that results in product damage, inaccurate data or the cancellation of a procedure.  
**Notes**  
**A Note** proves useful information in terms of a function or procedure.  
**Security Measures**

CAUTIONS

Federal (US) law restricts this device to sale by or on the orders of a physician.

**Notes:**  
This manual is written to protect the maximum settings. Therefore, your model may or may not have the indicators and functions described, depending on what you ordered.  
 This unit uses internal power and has an IEC.EN 60601-1 Type BF safety. The type of protection BF means that the connection between the unit and the user meets the electrical and dielectric leakage safety standards of IEC.EN 60601-1. **WARNING** and **CAUTION** messages must be followed. To avoid possible injury, follow all instructions when operating the unit.

CAUTIONS


- Only used by health practitioners on doctor's orders.
- Before the DP1 is determined for personal use, the user must receive proper instructions or training.
- Ultrasonic Pocket Fetal Doppler does not demonstrate preterm delivery or prevention of preterm birth.
- The Doppler unit is a tool for healthcare practitioners and should not be used in the setting of normal fetal detection. It is not intended for treatment.
- Placement of the ultrasound transducer on the abdomen is essential for obtaining the fetal heartbeat as opposed to the mother's heartbeat or sounds from the abdomen. Users must receive proper unit placement technique training either through acceptable Ob/Gyn training and individual circumstances or training from a physician and have been trained in unit placement.
- The unit is not explosion proof and cannot be used near flammable objects.
- Magnetic and electric fields can affect the unit's ability to properly. For this reason, ensure that all devices operating in the vicinity of the unit comply with EMC regulatory requirements. X-ray equipment and magnetic resonance imaging (MRI) devices can emit high levels of electromagnetic radiation.
- We recommend that ultrasound exposure be kept as low as possible. This is really considered good practice and should be observed at all times.
- Do not use the unit in conjunction with high-frequency surgical instruments and do not use the unit in the vicinity of magnetic resonance imaging (MRI) equipment which can emit high electromagnetic radiation.
- The unit is not protected against defibrillation.
- ELECTRIC SHOCK HAZARD-** Do not attempt to change the battery when your hands are wet.
- Do not connect any equipment or accessories that are not approved by the manufacturer or do not meet IEC 60601-1 standards to the unit. Use for unauthorized equipment or accessories has not been tested or supported and operation of the unit and safety are not guaranteed.
- The use of accessories outside of the manufacturer's specifications may result in increased electromagnetic emissions or reduced electromagnetic resistance of the unit.
- Units must not be used adjacent to or stacked on top of each other and if adjacent or stacked use is required, the units must be considered to ensure the operation to which they are to be used.
- Electro-medical equipment requires installation and is put into service in accordance with the EMC information provided in this manual.
- Radio communication equipment whether portable or cellular affects electrical equipment. See the **Recommended Distance section**.
- Do not repair or maintain any device or accessory that is being used with a patient.

CAUTIONS

- Provide service to professional technicians.
- Store the device in a clean environment and avoid vibration during storage.
- Do not sterilize Doppler with steam or gas.
- Electromagnetic Interference** - Ensure that the environment in which the device is operated is not adjacent to strong electromagnetic sources, such as radio transmitters, cellular telephones, etc.
- Prior to a Doppler examination, check for any visible damage to the unit and probe that could endanger the patient/operator or machine performance. If damage is found, replace it with a decent one.
- The following safety checks must be carried out every two years or as determined in the testing and inspection process by a qualified person who has adequate training, knowledge and practicality to carry out these tests.
  - Check equipment for mechanical and functional damage.
  - Check the security label.
  - Check equipment for mechanical and functional damage.
  - Check security labels that are relevant and easy to read.The leakage current should not exceed the limit. Data must be recorded in unit records. If the unit does not function properly or fails during testing, the unit must be repaired.
- Devices and accessories must be disposed of in accordance with regulations after their lifetime, or may be returned to a dealer or recycling plant or appropriate disposal according to hazardous waste. Do not dispose of the unit together with household waste.

**Introduction**  
**Purpose of Use / Indications for Use**  
The DP1 is a device intended for detecting the fetal heartbeat. The DP1 is intended to be operated by a professional and has been trained by a doctor for use by pregnant women in hospitals, clinics or at home.

- Feature**
- Detect and display FHR
  - Fetal heartbeat sound Heartbeat
  - Auto power off in ±2 minutes
  - Sound level setting button
  - Sound level setting
- Heart rate signal intensity icon  
Battery indicator  
Low battery warning



LCD Display & Touch Keys



Items	Description	
1	Heart rate	Heart rate indicator and flashes when it detects a heartbeat
2	Heart rate signal intensity indicator	This indicator is displayed on the left side of the screen and has three states: empty, half empty and full. Which respectively represent low, medium and high fetal heart signal intensity.
3	Heart rate indicator	Displays the fetal heart rate in the range from 50 bpm to 240 bpm. When the fetal heart rate is out of range,
	Volume indicator	The volume number is displayed in the center of the screen, the same area as the FHR Numeric. The volume display setting will appear for 0.5 seconds, then the display will return to the FHR numeric display. Sound volume ranges from 0~7.
4	Battery Indicator	The battery indicator is displayed on the right side of the screen. There are 5 battery levels, represented by 0~\$ panels. If the battery is empty then the icon will flash.
5	Volume increase touchpad	Touch and hold the button briefly to increase the volume.
6	Volume reduction touchpad	Touch and hold the button briefly to decrease the volume.
7	Touchpad On/Off	When doppler is off, touch and hold this button to turn on. When doppler is on, touch and hold this button to turn it off.

**Battery**  
The DP1 uses 2 AA lithium batteries. Battery specifications LR6, AA 1.5V  
**Notes:**  
You can use AA lithium batteries for the same specifications and can be purchased at the nearest store.

**Basic Operation**  
**Notes:**  
To ensure the doppler is working properly, please read this section and the **Precautions section** before operating; follow the steps when connecting all components.  
**Opening and Checking Packages**  
Open the Package, carefully remove the Doppler and accessories. Take care of the packaging for transportation or storage at a later date. Check the components with the Packing List.  
Check all External Damage.  
Check all cables and accessories.  
If there is a problem, contact you or the distributor directly.  
**Battery Installation.**  
a. Remove the screw with a plus screwdriver and remove the battery cover.  
b. Insert the battery into the battery compartment carefully. Make sure the positive and negative poles are installed correctly, for installation can be seen in the battery compartment.  
c. Install the battery cover and lock with screws.

**Battery removal/Installation**  
a. Remove the screw with a plus screwdriver and remove the battery cover  
b. Take the battery that has been used, you can immediately replace it with a new one, make sure the battery is installed correctly  
c. Install the battery cover and lock with screws

CAUTIONS

- Turn off Doppler before removing battery
- Replace alkaline batteries with batteries provided by the manufacturer or buy in stores according to the required specifications. See the product specifications section for detailed battery specifications.
- If the battery is not installed properly, the Doppler will not work properly or it will be damaged
- Do not disassemble the shorted battery
- Do not recharge the battery
- Do not throw the battery in fire or water
- Do not stick metal objects on the battery circuit
- Do not combine and use batteries with different types
- Do not solder the battery directly. If you need soldering and welding, consult our technician for a suitable method
- Don't overuse the battery
- To replace the battery, follow the product instructions for use
- Keep the battery away from children. If swallowed, contact a doctor immediately
- Store the battery in a dry and cool place. Do not store the battery at a temperature of 45°C and above or at a humidity of 75% and above
- Dispose of batteries according to regulations according to IEC61429 for standard disposal if necessary
- Remove and store the battery in a cool, dry place if the Doppler is not used for a long time
- A battery that has a lifetime. If the Doppler battery usage becomes shorter than normal, the battery life has expired. Replace the battery with a new battery with the same specifications as the manufacturer.

**Turn on**  
Press the On/Off button for a few moments while the Doppler is off and the Doppler will display a lit display before switching to a test display  
**turn off**  
Press the On/Off button for a while while Doppler is on and Doppler will be off. If Doppler is not used for 2 minutes, Doppler will turn off automatically  
**FHR detection**  
Before using Doppler to look for FHR, you should always make sure the Doppler is in good condition and if any damage will affect patient safety and the function of the device. If damage is found, discontinue use and replace with a new one.

**Procedure for detecting FHR**  
a. Patient lying face up  
b. Apply the gel on the surface of the ultrasonic Doppler transducer and turn on the Doppler  
c. Touch the patient's abdomen to determine the location of the fetus  
d. Apply Doppler on the patient's abdomen and rotate it around the fetus until a clear heart sound is obtained and the FHR number is stable

**Notes:**  
1. Don't mistake the mother's heart rate with the fetal heart rate  
2. Do not use gloves to touch the screen. If there is water and gel on the finger, please clean it first or touch will have an effect

The best way to find a heartbeat signal

- The easiest and fastest way: refers to the position of the fetal heart at the last position detected by the doctor as a reference and Doppler movement around that position slowly until the best FH signal is found
- The position of the fetal heart can change when the fetus moves in the uterus. You can confirm the position of the fetus in advance according to the uterine fundus (top of the uterus) at different weeks of pregnancy
  - At the end of the 12th week of gestation, the height of the uterine fundus is 2-3 fingers above the symphysis pubis (about 2-3 cm).
  - At the end of the 16th week of gestation, the height of the uterine fundus is midway through the symphysis and the center of the uterus
  - At the end of the 20th week of gestation, the height of the uterine fundus is 1 cm below the center
  - At the end of the 24th week of gestation, the height of the uterine fundus is 1 cm above the center of the uterus
  - At the end of the 28th week of gestation, the height of the uterine fundus is 3 cm above the center
  - At the end of the 32nd week of gestation, the height of the uterine fundus is between the xiphisternum and the center of the uterus
  - At 36 weeks of gestation, the height of the uterine fundus is 3 cm below the xiphisternum



The clarity and volume level of the fetal heartbeat is obtained from the fetal back. Fetal movement is usually the movement of the fetal limbs. So, if the movement of the fetus is often obtained from the right abdomen. The possibility of the fetal back is on the left and vice versa. You can find the fetal back based on the fetal movement  
If the birth is head down (cephalic), the fetal heart is between the right or left of the center  
**Steps to Find Fetal Heartbeat**  
Position the patient in a supine and relaxed position >> Confirm the fetal position by hand >> Apply a little gel to the Doppler >> Place the Doppler on the patient's abdomen and start looking for the fetal heart >> The fetal heart can be found when the Doppler sounds "Boom-Boom-Boom "

CAUTIONS

- The Doppler has an IP22 protection system that is water resistant. Do not drop the unit into the water
- Doppler is very subtle and sensitive. It is expected to be careful when carrying to avoid falling on the ground or other hard surfaces. All forms of damage caused by falling are not covered by the warranty
- Keep the gel away from children. If swallowed, contact a doctor

**Notes:**  
1. The best heart rate quality is only obtained when the Doppler is placed where the signal is best detected  
2. Do not place the Doppler near the placenta or the central bloodstream  
3. If the fetus is in a head-down (cephalic) position and the mother is in a supine position, the sound of the heartbeat is most pronounced in the middle below the navel. During the examination, prolonged supine position should be avoided, to avoid high blood pressure. Providing a bolster or pillow under the patient's head or feet can help reduce this risk  
4. If you do not get an optimal heart rate signal, it does not guarantee the accuracy of the reading. If the FHR (Fetal Heart Rate) reading does not match the sound of the heartbeat, the heartbeat sound from the reading will prevail  
5. When worn on a patient, the Doppler temperature may feel slightly warm (<2°C(35.6°F)) above room temperature. When not in use, the Doppler temperature may be slightly (<5°C(41°F)) above room temperature.

**After Use**  
1. Turn off Doppler  
2. Clean the remaining gel from the patient and probe with a soft cloth or tissue

**Care and Cleaning**  
**Care**  
Before use, check the completeness and safety of the user or Doppler function. In case of damage, contact the manufacturer for immediate service or replacement  
All Doppler checks, including function and safety checks, must be carried out by trained personnel, every check is carried out every 12 months and after every service. And the safety check should include leakage current test and insulation test. In addition to the above requirements, comply with local regulations regarding maintenance and measurement. The accuracy of the FHR reading is determined by Doppler and can be set by the user, if in doubt about the FHR reading, take the measurement by another method, such as a stethoscope or contact your local distributor or manufacturer for assistance. Doppler is easily damaged and must be handled with care. Clean the rest of the gel from the Doppler after each use. Regular use of the unit, can help extend the life of the Doppler  
Replace accessories such as batteries when they run out. If one of the accessories is damaged, read the Product information section for details and buy a new one

**Cleaning**  
Before cleaning, turn off Doppler. Keep the surface of the body clean and free from dust and dirt. Clean the body with a soft dry cloth. If possible clean with a cloth dampened in soapy water, ethanol (75%) or isopropanol (70%). Then wipe with a dry cloth.

CAUTIONS

- Do not use harsh solvents such as acetone
- Never use objects that can steal wood or steel brushes
- The Doppler is protected and has splash protection with an IP22 rating. Do not drop the Doppler directly into the water
- Do not give any solution after cleaning the surface

**Disinfection**  
On normal surfaces, Doppler does not need disinfection. If the surface is dirty, clean the unit and then disinfect the unit with a soft cloth and dampen it with ethanol (75%) or isopropanol (70%). Then wipe with a dry cloth

CAUTIONS

Do not dampen Doppler in a disinfectant

**Sterilization**  
Do not sterilize Doppler  
**Notes:**  
After cleaning or disinfection, check the Doppler function. If any problems are detected, contact the manufacturer or service before using again

Item Check	Method
Visual Check	Check Doppler for any damage
Function check	Check Doppler, whether the bias is turned on normally or not (See how to turn it on and off). When Doppler is turned on, check whether the display panel is lit normally like <b>LCD Display and Touch Key</b> : touch the ultrasound transducer head and check whether it can produce clear sound or not

Product Specification	
Product Information	
Product name	Ultrasonic Pocket Doppler
Model	DP1

Complied Standards	
IEC 60601-1:2005, EN 60601-1:2006/A1:2003, IEC 60601-1-2:2004<IEC 60601-2-37:2015, EIC 61266:1994	

Classification	
Anti Shock Type	AA Battery Power
Anti Shock Class	Type BF
Water Resistance Class	IP22, resistance to drops or splashes when the BATTERY COVER is tilted 15°
Resistance To Flammable Gases Class	The unit cannot be used near flammable gas
Work system	Continuously
EMC (Electromagnetic Compatibility)	CISPR 11 group 1 Class B




Physical Specification	
Size	Length x Width x Height: 48 mm(P) x 39 mm(L) x 147 mm(T)
Weight	90 g (excluded battery)
LCD	size 1.1" inch
	Display <ul style="list-style-type: none"><li>FHR</li><li>Battery Level</li><li>Signal Intensity</li></ul>
Gel	pH: 5.5~8.0 Acoustic Impedance: 1.5 x 10 <sup>6</sup> Pa.s/m ~1.7 x 10 <sup>6</sup> Pa.s/m (35°C/95°F)

Environmental conditions	
Working:	Temperature: +5°C-40°C Humidity: 15%-80% Atmosphere Pressure : 700hPa
Transport and storage	Temperature: -20°C-55°C Humidity: ≤95% Atmosphere Pressure : 700hPa

**Note:** the time it takes Doppler to use from the minimum temperature is 2 hours. The time required for Doppler to cool from the maximum storage temperature until it is ready for use is 2 hours

Product Specification	
FHR (fetal heart rate)	Range : 50bpm~240bpm Accuracy: ±3 bpm Note: FHR measurement results may not be accurate if Doppler measurement is out of range
FHR Resolution	1 bpm
Audio Output	Output: 2W, Noise: ≤45 dBa
Overall sensitivity	>110 dB
Auto Off	The unit will turn off automatically if the Doppler does not receive a signal or is not operated for 2 minutes
Bluetooth	Transmission distance (without obstruction) >5m (inside the room depending on the building structure and material
Ultrasound	Normal Frequency: 3 Mhz Working Frequency: 3 Mhz P_<1 MPa, Lob <10 mW/cm <sup>2</sup> Ispta < 100 mW/cm <sup>2</sup> , Isata < 10mW/cm <sup>2</sup> Isppa.3 < 190 mW/cm <sup>2</sup> , Ispta.3 <94mW/cm <sup>2</sup>



7		Serial Number	16	Rx Only	US law restricts this product from being sold by a doctor's recommendation
8		Manufacture Date	17	FCC ID: SMQSD1ME AND	Federal Communications Commission: FCC ID: SMQSD1ME AND
9		Manufacturer			

CONTROLLED COPY